

Kinect-Based Human Motion Analysis

M.Sc. David Enrique Veloz Renteria

Contents

- ▶ Experimental Setup
- ▶ Flowchart
- ▶ UI
- ▶ Software Features
- ▶ Test Setup
- ▶ Standing Up and Sitting Down
- ▶ Posture with Feet Together and Eyes Open and Closed
- ▶ Walking in Place
- ▶ Results


Experimental Setup

► Materials:



- Tripod
- Kinect V2
- Computer running Windows 11
- MATLAB



Kinect on a tripod




Flowchart



UI



Main interface





SAS interface

Software Features


- ▶ **Frei:** The software detects a person and tracks the 25 joints identified by Kinect
- ▶ **Standing and Sitting (SAS):** detects whether the person is sitting or standing
- ▶ **Posture with feet together and eyes open and closed (POCO):** detects when the person brings their feet together or apart; the eye status of the person monitoring the test must be confirmed
- ▶ **Walking in place (SIP):** detects an upright person to begin data recording

Test Setup

- ▶ 5 healthy subjects: 2 women and 3 men (students and staff of TU Ilmenau)
- ▶ 3 SAS
 - ▶ Standard chair height
 - ▶ Comfortable chair height
 - ▶ Max. chair height
- ▶ POCO
- ▶ SIP
- ▶ Verify the validity of Kinect for analyzing human movements



Standing Up and Sitting Down



Test time [s]	55
Number of repetitions:	5
Average displacement [m]:	0.3662
Seat height [cm]	46
↑ - Time [s]	0.9533
↑ Range of motion: medial-lateral [cm]	4.3742
↑ Anterior-posterior range of motion [cm]	43.1813
↑ Range of motion of the right hand [cm]	27.648
↑ Range of motion of the left hand [cm]	27.6949
↓ - Time [s]	0.8667
↓ Range of motion: medial-lateral [cm]	3.3855
↓ Range of motion: Anterior-posterior [cm]	40.5402


Posture with Feet Together and Eyes Open and Closed




Test time [s]	120
Pitch deflection range [°]	1.0586
Roll deflection range [°]	0.4296
3D deflection range [°]	3.0946
Average pitch fluctuation rate [°/s]	0.0529
Average roll fluctuation rate [°/s]	0.0215
Average 3D fluctuation rate [°/s]	0.1547



Pitch deflection range [°]	0.3751
Roll deflection range [°]	0.4774
3D deflection range [°]	-4.4292
Average pitch fluctuation rate [°/s]	0.0188
Average roll fluctuation rate [°/s]	0.0239
Average 3D fluctuation rate [°/s]	-0.2215



Walking in Place



Test time [s]	125
Range of motion, left knee [cm]	19.81
Range of motion, right knee [cm]	18.63
Cadence left [steps/min]	43.4694
Cadence right [steps/min]	43.299

Results

	Gakivax		Kinect (Paper)		Vicon (Paper)	
	Mittelwert	SA	Mittelwert	SA	Mittelwert	SA
Aufstehen und Hinsetzen						
Durchschnittliche Verschiebung [m]:	0.32	0.06	-	-	-	-
↑ - Time [s]	1.28	0.92	1.31	0.25	1.38	0.25
↑ Range of motion: medial-lateral [cm]	2.64	2.40	2.42	1.42	2.15	1.16
↑ Anterior-posterior range of motion [cm]	34.45	6.30	37.02	6.51	37.22	6.41
↑ Range of motion of the right hand [cm]	25.27	9.69	31.74	6.28	32.99	6.34
↑ Range of motion of the left hand [cm]	23.53	7.95	34.04	6.96	34.67	7.07
↓ - Time [s]	1.72	1.45	1.53	0.23	1.59	0.23
↓ Range of motion: medial-lateral [cm]	1.79	2.42	2.87	1.45	2.37	1.26
↓ Range of motion: Anterior-posterior [cm]	30.55	11.33	39.67	6.27	39.29	5.97
Haltung mit geschlossenen Füßen und offenen und geschlossenen Augen						
EO Pitch deflection range [°]	0.75	0.30	1.67	0.76	1.49	0.63
EO Roll deflection range [°]	0.76	0.28	1.54	0.52	1.35	0.4
EO 3D deflection range [°]	0.04	0.01	0.17	0.06	0.18	0.06
EO Average pitch fluctuation rate [°/s]	0.04	0.01	0.21	0.07	0.2	0.06
EC Pitch deflection range [°]	0.58	0.28	1.67	0.76	1.49	0.63
EC Roll deflection range [°]	0.81	0.31	1.54	0.52	1.35	0.4
EC 3D deflection range [°]	0.03	0.01	0.17	0.06	0.18	0.06
EC Average pitch fluctuation rate [°/s]	0.04	0.02	0.21	0.07	0.2	0.06
Gehen auf der Stelle						
Range of motion, left knee [cm]	14.34	6.34	4.37	1.06	9.05	1.7
Range of motion, right knee [cm]	14.09	6.47	4.38	1.02	9.03	1.89
Cadence left [steps/min]	49.02	10.73	47.08	7.31	47	7.21
Cadence right [steps/min]	46.56	10.08	47.63	7.18	47.58	7.15




Results (2)

- ▶ The algorithms used to calculate the individual parameters are not explained in the researched papers
- ▶ Development of custom algorithms
- ▶ Pseudocode for “marching in place”:

Start

```
|| Read knee signals (t, x, y, z)
|| Set thresholds for both knees
|| Search for maximum values in the original knee signals
|| Fit a second-degree polynomial trend curve to the found maximum values
|| Use the trend curve to normalize the data and set new thresholds
|| Search for maximum values in the normalized signals of both knees
|| Calculate the parameters
|| Print results
```

End



Results (3)

- ▶ Summary of the results in a report
- ▶ Report:
 - ▶ key charts for each test
 - ▶ calculated numerical values
 - ▶ manually done

